

Product Datasheet

Canine TWEAK (TNFSF12) Recombinant Protein (orb1820612)

Description	The Canine TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	TWEAK
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NASKGRKTRA RRAIAAHYEV HPQPGQDGGQ AGVDGTVSGW EEAKINSSNP LRYDRQSGEF IVTRAGLYYL YCQVHFDEGK AVYLKLDLLV DDALALRCLE EFSATAASSL GPQLRLCQVS GLLPLRPGSS LRIRTLPAWAH LKAAPFLTYF GLFQVN (156)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	100568280
Expiration Date	6 months from date of receipt.