

## Product Datasheet

# Mouse IL-1 Receptor Antagonist Recombinant Protein (orb1820616)

<b>Description</b>	The Mouse IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally ...
<b>Target</b>	IL-1RA
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SGKRPCMQA FRIWDTNQKT FYLRNNQLIA GYLQGPNIKL EEKIDMVPID LHSVFLGIHG GKLCLSCAKS GDDIKLQLEE VNITDLSKNK EEDKRFTFIR SEKGPTTSFE SAACPGWFLC TTLEADRPVS LTNTPEEPLI VTKFYFQEDQ (150)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.1 kDa
<b>Entrez</b>	<a href="#">16181</a>
<b>Expiration Date</b>	6 months from date of receipt.