

Product Datasheet

Swine BAFF (TNFSF13B) Recombinant Protein (orb1820617)

Description	The Swine BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAQDAEETVT QDCLQLIADS DMPTIRKGAY TFVPWLLSFK RGRALEEKEN KIVVKETGYF FIYGQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISRD GDGTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	100038026
Expiration Date	6 months from date of receipt.