

## Product Datasheet

# Swine TWEAK (TNFSF12) Recombinant Protein (orb1820624)

<b>Description</b>	The Swine TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	TWEAK
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SAPKGRKPRA RRVIAAHYEV HPRPGQDGAQ AGVDGTVSGW EEAKINSSNP LRYDRQSGQF LVTRAGLYYL YCQVHFDEGK AVYLKLDLLV DDTLALRCLE EFSATAASSL GPQLRLCQVS GLLPLRPGSS LRIRLPWAH LKAAPFLTYF GLFQVH (156)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.3 kDa
<b>Entrez</b>	<a href="#">100127163</a>
<b>Expiration Date</b>	6 months from date of receipt.