

Product Datasheet

Swine TWEAK (TNFSF12) (orb1820626)

Description	The Swine TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	TWEAK
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SAPKGRKPRA RRVIAAHYEV HPRPGDGAQ AGVDGTVSGW EEAKINSSNP LRYDRQSGQF LVTRAGLYYL YCQVHFDEGK AVYLKLDLLV DDTLALRCLE EFSATAASSL GPQLRLCQVS GLLPLRPGSS LRIRTLPWAH LKAAPFLTYF GLFQVH (156)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	100127163
Expiration Date	6 months from date of receipt.