

Product Datasheet

Canine IL-1 alpha Recombinant Protein (orb1820631)

Description	The Canine IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVAYNFHNNE KYNIRIIS QFILNDNLNQ SIVRQTGGNY LMTAALQNLNLD DAVKFDMGAY TSEDSKLPVT LRISKTRLFV SAQNEDEPVL LKEMPETPKT IRDETLLFF WERHGSKHYF KSVAQPKLFI ATQERKLVHM ARGQPSITDF RLLETQP (157)
Purity	98%
Source	Yeast
MW	18.2 kDa
Entrez	403782
Expiration Date	6 months from date of receipt.