

## Product Datasheet

# Canine BAFF (TNFSF13B) Recombinant Protein (orb1820637)

<b>Description</b>	The Canine BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	BAFF
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AIQGPEETVT QDCLQLIADS DTPTIRKGAY TFVPWLLSFK RGRALEEKEN KILVKEAGYF FIYSQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISRD GDGTFFGALK LL (152)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.1 kDa
<b>Entrez</b>	<a href="#">485545</a>
<b>Expiration Date</b>	6 months from date of receipt.