

Product Datasheet

Canine BAFF (TNFSF13B) Recombinant Protein (orb1820638)

Description	The Canine BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AIQGPEETVT QDCLQLIADS DTPTIRKGAY TFVPWLLSFK RGRALEEKEN KILVKEAGYF FIYSQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISRD GDGTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	485545
Expiration Date	6 months from date of receipt.