

Product Datasheet

Bovine BAFF (TNFSF13B) Recombinant Protein (orb1820640)

Description	The Bovine/Ovine/Caprine BAFF yeast-derived recombinant protein is not tagged and is naturally endot...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ALQEAETVT QDCLQLIADS DTPTIRKGAY TFVPWLLSFK RGRALEEKEN KILVKETGYF FIYGQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISRD GDGTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	504507
Expiration Date	6 months from date of receipt.