

Product Datasheet

Canine IFN alpha Recombinant Protein (orb1820641)

Description	The Canine IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CHLPDTHGLR NWRVLTLLGQ MRRLSAGSCD HYTNDFAFPK ELFDGQRLQE AQALSVVHVM TQKVFHLFCP DTSSAPWNMT LLEELCSGLS EQLDDLEACP LQEAGLAETP LMHEDSTLRT YFQRISLYLQ DRNHSPCAWE MVRAEIGRSF FSSTILQERI RRRK (164)
Purity	98%
Source	Yeast
MW	19.0 kDa
Entrez	100682977
Expiration Date	6 months from date of receipt.