

Product Datasheet

Rat IL-17A Recombinant Protein (orb1820646)

Description	The Rat IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVLIPQSSVC PNAEANNFLQ NVKVNKVLN SLSSKASSRR PSDYLNKSTS PWTLNRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHHMNS VLIQQEILVL KREPEKCPFT FRVEKMLVGV GCTCVSSIVR HAS (133)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	301289
Expiration Date	6 months from date of receipt.