

Product Datasheet

Rabbit IL-17A Recombinant Protein (orb1820648)

Description	The Rabbit IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIAMPRNPGC PNAEDKNFPQ NVKVSLNILN KSVNSRRPSD YYNRSTSPWT LHRNEDHERY PSVIWEAKCR HLGCVNAEGN VDHHMNSVPI QQEILVLRRE SQHCPHSFRL EKMLVAVGCT CVTPIIHHMA (130) [Arg57His]
Purity	98%
Source	Yeast
MW	14.9 kDa
Entrez	100136642
Expiration Date	6 months from date of receipt.