

Product Datasheet

Dolphin IL-1 Receptor Antagonist Recombinant Protein (orb1820655)

Description	The Dolphin IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
Target	IL-1ra
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	YPLGKRPCEM QAFRIWDVNO KTFYLRNNQL VAGYLQGPNT KLEEKIDVVP IEPHAMFLGI HGGKLCLACV KSGDEIKLGL EPVNITDLNS SKEEDKRFAF IRSDSGPTTS FESAACPGWF LCTALETDPQ VGLTNTPDQA VQVTKFYFQQ DQ (152)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	101328434
Expiration Date	6 months from date of receipt.