

## Product Datasheet

### Human IL-19 Recombinant Protein (orb1820657)

<b>Description</b>	The Human IL-19 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	IL-19
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LRRCLISTDM HHIEESFQEI KRAIQAKDTF PNVTILSTLE TLQIIKPLDV CCVTKNLLAF YVDRVFKDHQ EPNPKILRKI SSIANSFLYM KQTLRQCQEQ RQCHCRQEAT NATRVIHDNY DQLEVHAAAI KSLGELDVFL AWINKNHEVM FSA (153)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.8 kDa
<b>Entrez</b>	<b>29949</b>
<b>Expiration Date</b>	6 months from date of receipt.