

Product Datasheet

Mouse TNFSF11 (RANKL) Recombinant Protein (orb1820660)

Description	The Mouse TNFSF11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	TNFSF11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FSGAPAMMEG SWLDVAQRGK PEAQPF AHLT INAASIPSGS HKVTLSSWYH DRGWAKISNM TLSNGKLRVN QDGFYYLYAN ICFRHHETSG SVPTDYLQLM VYVVKTSIKI PSSHNLMKGG STKNWSGNSE FHFYSINVGG FFKL RAGEEI SIQVSNP SLL DPDQDATYFG AFKVQDID (178)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	21943
Expiration Date	6 months from date of receipt.