

Product Datasheet

Rabbit CCL4 Recombinant Protein (orb1820663)

Description	The Rabbit CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTA CCFSYTLRKL PRHFVIDYFE TSSLCSQPAV VFQTKKGRQV CANPSES WVQ EYVDDLELN (69) [Thr32Ser]
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	100008984
Expiration Date	6 months from date of receipt.