

Product Datasheet

Mouse IL-4 Recombinant Protein (orb1820668)

Description	The Mouse IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HIHGCDKNHL REIIGILNEV TGEGTPCTEM DVPNVLATK NTTESELVCR ASKVLRI FYL KHGKTPCLKK NSSVLMELQR LFRAFRC LDS SISCTMNESK STSLKDFLES LKSIMQMDYS (120)
Purity	98%
Source	Yeast
MW	13.6 kDa
Entrez	16189
Expiration Date	6 months from date of receipt.