

Product Datasheet

Mouse IL-4 Recombinant Protein (orb1820669)

| | |
|-------------------------|--|
| Description | The Mouse IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | IL-4 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | HIHGCDKNHL REIIGILNEV TGEGTPCTEM DVPNVLATK NTTESELVCR ASKVLRI FYL KHGKTPCLKK NSSVLMELQR LFRAFRC LDS SISCTMNESK STSLKDFLES LKSIMQMDYS (120) |
| Purity | 98% |
| Source | Yeast |
| MW | 13.6 kDa |
| Entrez | 16189 |
| Expiration Date | 6 months from date of receipt. |