

Product Datasheet

Rabbit TNF alpha Recombinant Protein (orb1820670)

Description	The Rabbit TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SASRALSDKP LAHVVANPQV EGQLQWLSQR ANALLANGMK LTDNQLVVPA DGLYLIYSQV LFSGQGCRSY VLLTHTVSRF AVSYPNKVNL LSAIKSPCHR ETPEEAEPMA WYEPIYLGGV FQLEKGDRLS TEVNQPEYLD LAESGQVYFG IIAL (154)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	100009088
Expiration Date	6 months from date of receipt.