

Product Datasheet

Canine IL-1 Receptor Antagonist Recombinant Protein (orb1820674)

Description	The Canine IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-1ra
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LGKRPCRMQAFRIWDVNQKTFYLRNNQLVAGYLQGSNTKLEEKLDVVPVEPHA VFLGIHGGKLCACVKSGDETRLQLEAVNITDLSKNKDQDKRFTFILSDSGPTTSF ESAACPGWFLCTALEADRPVSLTNRPEEAMMVTKFYFQKE (149)
Purity	98%
Source	Yeast
MW	16.8 kDa
Entrez	403660
Expiration Date	6 months from date of receipt.