

## Product Datasheet

# Swine IL-1 Receptor Antagonist Recombinant Protein (orb1820680)

<b>Description</b>	The Swine IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally ...
<b>Target</b>	IL-1ra
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HPLGKRPCRQ QAFRIWDVNO KTFYLRNNQL VAGYLQGPNT KLEEKIDVVP VEPHFVFLGI HGGKLCCLSCV KSGDEMKLQL EPVNITDLRK NSEQDKRFTF IRSDSGPTTS FESAACPGWF LCTALEADQP VGLTNTPKAA VKVTKFYFQQ DQ (152) [AspAla81/82GluPro}
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.2 kDa
<b>Entrez</b>	<a href="#">397499</a>
<b>Expiration Date</b>	6 months from date of receipt.