

Product Datasheet

Swine IL-1 Receptor Antagonist Recombinant Protein (orb1820681)

Description	The Swine IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally ...
Target	IL-1ra
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HPLGKRPCR M QAFRIWDV NQ KTFYLRNNQL VAGYLQGPNT KLEEKIDVVP VEPHFVFLGI HGGKLCCLSCV KSGDEMKLQL EPVNITDLRK NSEQDKRFTF IRSDSGPTTS FESAACPGWF LCTALEADQP VGLTNTPKAA VKVTKFYFQQ DQ (152) [AspAla81/82GluPro}
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	397499
Expiration Date	6 months from date of receipt.