

Product Datasheet

Dolphin IL-4 Recombinant Protein (orb1820683)

Description	The Dolphin IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKCDVTLQEI IKTLNILTAK KNLCMELPVE DVFATTKNTT EKETFCRAGT VLRHIYRHHK CFNKPLSGLH RNLSSMANMT CSVNEAKKST LKDFLERLKM IMKEKYSKC (109) [Gln64Lys}
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	101338293
Expiration Date	6 months from date of receipt.