

Product Datasheet

Mouse BAFF Recombinant Protein (orb1820692)

Description	The Mouse BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AFQGPEETEQ DVDLSAPPAP CLPGCRHSQH DDNGMNLRLNI IQDCLQLIAD SDTPTIRKGT YTFVPWLLSF KRGNALEEKE NKIVVRQTGY FFIYSQVLYT DPIFAMGHVI QRKKVHVFGD ELSLVTLFRC IQNMPKTLPN NSCYSAGIAR LEEGDEIQLA IPRENAQISR NGDDTFFGAL KLL (183)
Purity	98%
Source	Yeast
MW	20.6 kDa
Entrez	24099
Expiration Date	6 months from date of receipt.