

Product Datasheet

Mouse IFN gamma Recombinant Protein (orb1820703)

Description	The Mouse IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HGTVIESLESLNYYFNSSGIDVEEKSLFLDIWRNWQKDGMKILQSQIISFYLRLE EVLKDNQAISNNISVIESHLITFFSNSKAKKDAFMSIAKFEVNNPQVQRQAFNELI RVVHQLLPESLRKRKRRC (133)
Purity	98%
Source	Yeast
MW	15.5 kDa
Entrez	15978
Expiration Date	6 months from date of receipt.