

Product Datasheet

Mouse IL-6 Recombinant Protein (orb1820705)

Description	The Mouse IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTSQVRRGDFTE D TTPNRPVYTTSQVGG LITHVLWEIVEMRKELCNGNSDCMN NDDALAENNLK LPEIQRNDG CYQTGYNQEICLLKISSGLLEYHSYLEYMKNNLKD NKKDKARVLQRDTETLIHIFNQEVKDLHKIVLPTPISNALLTDKLESQKEWLRTKTI QFILKSLEEF LKVTLRSTRQT (187)
Purity	98%
Source	Yeast
MW	21.7 kDa
Entrez	16193
Expiration Date	6 months from date of receipt.