

Product Datasheet

Mouse TNF alpha Recombinant Protein (orb1820708)

Description	The Mouse TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRSSSQNSSDKPVAHVVANHQVEEQLEWLSQRANALLANGMDLKDNLVPPA DGLYLVYSQVLFKGGQCPDYVLLTHTVSRFAISYQEKVNLLSAVKSPCKDTPEG AELKPWYEPIYLGGVFQLEKGDQLSAEVNLPKYLDFAESGQVYFGVIAL (156)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	21926
Expiration Date	6 months from date of receipt.