

Product Datasheet

Bovine CCL4 Recombinant Protein (orb1820730)

Description	The Bovine CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTACCFSYTLRKIPRNFVNDYFETSSLCSQPAVVFQTKKGRQVCANPS EPWVQEYVDDLELN (69)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	414347
Expiration Date	6 months from date of receipt.