

Product Datasheet

Mouse APRIL Recombinant Protein (orb1820732)

| | |
|-------------------------|--|
| Description | The Mouse APRIL yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | APRIL |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | AVLTQKHKKKHSLVHLVNVNITSKADSDVTEVMWQPVLRRGRGLEAQQGDIVRV WDTGIYLLYSQVLFHDVTFMTMGQVVSREGQGRRETLFRCIRSMPSDPDRAYNSC YSAGVFHLHQGDIITVKIPRANAKLSLSPHGTFLGFVKL (146) |
| Purity | 98% |
| Source | Yeast |
| MW | 16.4 kDa |
| Entrez | 69583 |
| Expiration Date | 6 months from date of receipt. |