

Product Datasheet

Human IL-8 Recombinant Protein (orb1820740)

Description	The Human IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EGAVLPRSAK ELRCQCIKTY SKPFHPKFIK ELRVIESGPH CANTEIIVKL SDGRELCLDP KENWVQRVVE KFLKRAENS (79)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	3576
Expiration Date	6 months from date of receipt.