

Product Datasheet

Human IL-4 Recombinant Protein (orb1820743)

Description	The Human IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKCDITLQEI IKTLNSLTEQ KTLCTELTVT DIFAASKNTT EKETFCRAAT VLRQFYSHHE KDTRCLGATA QQFHRHKQLI RFLKRLDRNL WGLAGLNSCP VKEANQSTLE NFLERLKTIM REKYSKCSS (129)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	3565
Expiration Date	6 months from date of receipt.