

## Product Datasheet

### Mouse IL-17A Recombinant Protein (orb1820746)

<b>Description</b>	The Mouse IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	IL-17A
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AAIIPQSSAC PNTEAKDFLQ NVKVNLKVFN SLGAKVSSRR PSDYLNRRST PWTLHRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHMMNS VLIQQEILVL KREPESCPFT FRVEKMLVGV GCTCVASIVR QAA (133)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	15.0 kDa
<b>Entrez</b>	<b>16171</b>
<b>Expiration Date</b>	6 months from date of receipt.