

## Product Datasheet

# Bovine IL-1 Receptor Antagonist Recombinant Protein (orb1820748)

<b>Description</b>	The Bovine IL-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	IL-1ra
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HPLGKRRCEM QAFRIWDVNO KIFYLRNNQL VAGYLQGPNT KLEEKIDVVP IEPHTMFLGI HGGKLCLACV KSGDEIKLKL EAVNITDLNQ NREQDKRFAF IRFDNGPTTS FESAACPGWF LCTSLEADQP VGLTNMPTEA LKVTKFYFQQ D (151)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.2 kDa
<b>Entrez</b>	<a href="#">281860</a>
<b>Expiration Date</b>	6 months from date of receipt.