

Product Datasheet

Guinea Pig CCL5 Recombinant Protein (orb1820753)

| | |
|-------------------------|---|
| Description | The Guinea Pig CCL5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free.... |
| Target | CCL5 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SPYASDTTPC CFAYISRALP RTHIKEYFYT SSKCSNLAVV FVTRKNRQVC ANPEKKWVRE YINSLEMS (68) |
| Purity | 98% |
| Source | Yeast |
| MW | 7.9 kDa |
| Entrez | 100135495 |
| Expiration Date | 6 months from date of receipt. |