

## Product Datasheet

# Dolphin IL-1 beta Recombinant Protein (orb1820754)

<b>Description</b>	The Dolphin IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-1 beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AAVQSLTCKL QDRDQKSLVL ASPCVLKALH LLARDMNREV VFCMSFVQGD ESNDKIPVAL GLKEKNLYLS CVMKGDRPIL QLEEVDPKTY PKWKMEKRFV FNKTEIKNSV EFESALYPNW YISTSQAEEK PIFLGRSKGG HDITDFTMEI ISP (153)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.5 kDa
<b>Entrez</b>	not available
<b>Expiration Date</b>	6 months from date of receipt.