

## Product Datasheet

# Rabbit IL-1 alpha Recombinant Protein (orb1820755)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Rabbit IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...  |
| <b>Target</b>           | IL-1 alpha  |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | SVPYTFQRNM RYKYLRIKQ<br>EFTLNDALNQ SLVRDTSQY<br>LRAAPLQNLG DAVKFDMGVY<br>MTSEDSILPV TLRISQTPLF VSAQNEDEPV<br>LLKEMPETPR IITDSEDIL FFWETQGNKN<br>YFKSAANPQL FIATKPEHLV<br>HMARGLPSMT DFQIS (155) |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 17.8 kDa  |
| <b>Entrez</b>           | <a href="#">100009250</a>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |