

Product Datasheet

Dolphin TNF alpha Recombinant Protein (orb1820769)

Description	The Dolphin TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRSSSKTSSN KPVAHVVANL STQGQLRWLN TYANTLLANS VKLEDNQLVV PTDGLYLIYS QVLFRGQGCP STHLFLTHTI SRIAVSYPTK VNLLSAIKSP CQRETPEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPD YLDAESGQV YFGIIAL (157) [Ser89Thr]
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	not available
Expiration Date	6 months from date of receipt.