

Product Datasheet

Swine IL-4 Recombinant Protein (orb1820774)

Description	The Swine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKCDITLQEI IKTLNILTAR KNSCMELPVT DVFAAPENTT EKETFCEAST VLRHIYRHHT CMKSLLSGLD RNLSSMANMT CSVHEAKKST LKDFLERLKT IMKEKYSKC (109)
Purity	98%
Source	Yeast
MW	12.5 kDa
Entrez	397225
Expiration Date	6 months from date of receipt.