

Product Datasheet

Feline IL-4 Recombinant Protein (orb1820783)

| | |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Description | The Feline IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | IL-4 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | QNFNNTLKEI IKTLNILTAR NDSCMELTVM DVLAAPKNTS DKEIFCRATT VLRQIYTHHN CSTKFLKGLD RNLSSMANRT CSVNEVKKCT LKDFLERLKA IMQKKYSKH (109) |
| Purity | 98% |
| Source | Yeast |
| MW | 12.6 kDa |
| Entrez | 751514 |
| Expiration Date | 6 months from date of receipt. |