

Product Datasheet

Feline IL-1 alpha Recombinant Protein (orb1820785)

Description	The Feline IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVAPNFYSSEKYNYQKIIKSQFILNDNLSQSVIRKAGGKYLAAAALQNLDDAVKFD MGAYTSKEDSKLPVTLRISKTRLFVSAQNEDEPVLLKEMPETPKTIRDETNLLFFW ERHGSKNYFKSVAHPKLFIAATQEEQLVHMARGLPSVTDFQILETQS (158)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	493944
Expiration Date	6 months from date of receipt.