

Product Datasheet

Swine IL-17F Recombinant Protein (orb1820788)

Description	The Swine IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-17F
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARKILKPGDP ALCPPLEDHS VRVDIRILSQ NQGVPLSHDF QNRSSSPWDY NITRDPHRFP SEIAEARCRH TGCINAQQQE DSSMNSVPIQ QEFLVLRREP QGCSRSFRLE KVRVTVGCTC VTPMVRYVRG LRAAGQLG (138)
Purity	98%
Source	Yeast
MW	15.5 kDa
Entrez	100037933
Expiration Date	6 months from date of receipt.