

Product Datasheet

Chicken IL-12 p40 Recombinant Protein (orb1820793)

Description	The Chicken IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally endotox...
Target	IL-12 p40
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AQWKLRENVY VIESEWNDET PAKKVKLTCD TSDEALPVYW KKGTELKGTG KTLTTEVKEF PDAGNYTCLS AKTHEIISYS FFLITKVDSN GQMIRSILKS YKEPSKTLK CEAKNYSGIF TCSWMTENES PSVKFTIRSL KGSQGDVTCS SPVARTDKSV TEYTAQCQKE NYCPFAEEHQ PTEMFLEVID EVEYENYTSS FFIRDIKPD PPQCQYASTN GTVTWYTPKT WSTPKSYFPL TFRVKVESTK KYKSKVYDAD EQSIQIPKTG PKDKISVQAR DRYYNSSWSE WSTLCR (296)
Purity	98%
Source	Yeast
MW	34.1 kDa
Entrez	404671
Expiration Date	6 months from date of receipt.