

Product Datasheet

Human IFN gamma Recombinant Protein (orb1820796)

Description	The Human IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QDPYVKEAEN LKKYFNAGHS DVADNGTLFL GILKNWKEES DRKIMQSQIV SFYFKLFKNF KDDQSIQKSV ETIKEDMNVK FFNSNKKKRD DFEKLTNYSV TDLNVQRKAI HELIQVMAEL SPAAKTGKRK RSQMLFRGRR ASQ (143)
Purity	98%
Source	Yeast
MW	16.8 kDa
Entrez	3458
Expiration Date	6 months from date of receipt.