

Product Datasheet

Canine IFN gamma Recombinant Protein (orb1820800)

| | |
|-------------------------|--|
| Description | The Canine IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | IFN gamma |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | QAMFFKEIEN LKEYFNASNP DVSDGGSLFV DILKKWREES DKTIIQSQIV SFYLKLFDFN KDNQIIQRSM DTIKEDMLGK FLNSSTSKRE DFLKLIQIPV NDLQVQRKAI NELIKVMNDL SPRSNLRKRK RSQNLFRGRR ASK (143) |
| Purity | 98% |
| Source | Yeast |
| MW | 16.7 kDa |
| Entrez | 403801 |
| Expiration Date | 6 months from date of receipt. |