

Product Datasheet

Rabbit IL-4 Recombinant Protein (orb1820802)

Description	The Rabbit IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GRRGDILPE VIKTLNILTE RKTPTKLM I ADALAVPKNT TEREAVCRAA TALRQFYLHH KVSWC FKEHG ELGDLRLLRG LDRNLCSMAK LSNCPGKEAR QTTLEDFLDR LKTAMQEKYS KRQS (124)
Purity	98%
Source	Yeast
MW	14.2 kDa
Entrez	100302454
Expiration Date	6 months from date of receipt.