

Product Datasheet

Canine TNF alpha Recombinant Protein (orb1820811)

Description	The Canine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VKSSSRTPSD KVAHVVANP EAEGQLQWLS RRANALLANG VELTDNQLIV PSDGLYLIYS QVLFKGQGCP STHVLLTHTI SRFAVSQTK VNLLSAIKSP CQRETPEGTE AKPWYEPIYL GGVFQLEKGD RLSAEINLPN YLDFAESGQV YFGIIAL (157)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	403922
Expiration Date	6 months from date of receipt.