

Product Datasheet

Canine IL-21 Recombinant Protein (orb1820814)

Description	The Canine IL-21 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-21
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKSSFQEQLL LLIRMRQLID IVDQLKNYVN DLDPESLPAP EDVKRHCERS AFSCFQKVQL KAANTGGNEQ IINVLTKQLK RKLPTNAGR RQKHPACPS CDSYEKAPPK EFLERLKS LI QKMIHQHLS (129)
Purity	98%
Source	Yeast
MW	14.9 kDa
Entrez	442935
Expiration Date	6 months from date of receipt.