

Product Datasheet

Rabbit IL-8 Recombinant Protein (orb1820827)

| | |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | The Rabbit IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | IL-8 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | AVRSLHCRLQ DAQQKSLVLS GTYELKALHL NAENLNQQVW FSMSFVQGEE SNDKIPVALG LRGKNLYLSC VMKDDKPTLQ LESVDPNRYK KKKMEKRFVF NKIEIKDKLE FESAQFPNWY ISTSQTEYMP VFLGNNSGGQ DLIDFSMEFV SS (152) |
| Purity | 98% |
| Source | Yeast |
| MW | 9.1 kDa |
| Entrez | 100009129 |
| Expiration Date | 6 months from date of receipt. |