

Product Datasheet

Rabbit IL-8 Recombinant Protein (orb1820828)

Description	The Rabbit IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVRSLHCRLQ DAQQKSLVLS GTYELKALHL NAENLNQQVW FSMSFVQGEE SNDKIPVALG LRGKNLYLSC VMKDDKPTLQ LESVDPNRYK KKKMEKRFVF NKIEIKDKLE FESAQFPNWY ISTSQTEYMP VFLGNNSSGGQ DLIDFSMEFV SS (152)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	100009129
Expiration Date	6 months from date of receipt.