

Product Datasheet

Rabbit IL-2 Recombinant Protein (orb1820833)

Description	The Rabbit IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTKETQEQLDQLLLDLQVLLKGVNDYKNSKLSRMLTFKFMPPKVKTELKH LQCLEELKPLEEVLNLAQGKNSHGGNTRESISNINVTVLKLGSETFMCEYDET VTIVEFLNRWITFCQSIISASSS (133)
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	100302458
Expiration Date	6 months from date of receipt.